

Labcorp Label Printing

07/24/2024 6:30 pm EDT

Labcorp Specimen Labels in DrChrono

In order to print Labcorp specimen labels, you must use a standard office printer. Labcorp labels do not print using a DYMO printer as Quest labels do.

Please use **Labcorp 8-Label paper** when printing Labcorp specimen labels. You may request this paper from your Labcorp sales representative. The Labcorp 8-Label paper is an 8.5x11" sheet of paper on which the order requisition prints at the top and 8 labels print at the bottom in the allocated spaces with peel-off labels.
